



Handicare 1000

A Minimum swivel radius	670mm
– with arm spacers	693mm
– with extended arm spacers	710mm
B Minimum wall to seat back	75mm†
C Seat depth	360mm
– with extended seat	425mm
D Minimum folded width to footplate	385mm†
E Minimum open width to edge of footplate	650mm†
F Armrest width – external	550mm
– with arm spacers	590mm
– with extended arm spacers	637mm
G Armrest width – internal	460mm
– with arm spacers	516mm
– with extended arm spacers	562mm
H Seat back height	453mm
I Footplate to seat height	465-540mm
– standard seat option	677-877mm
– perch seat option^	
J Minimum footplate height	85mm
Bottom parking space required	550mm (56°) min to 675mm at shallow angle (27°)††
Minimum track intrusion into staircase	115mm†
Minimum track overhang at top of stairs	86mm
Minimum track intrusion at bottom of stairs	400mm
Minimum track length with hinge	2300mm
Maximum weight limit	
– standard version	137kg/21.5 stone
– heavy duty version (XL)*	160kg/25 stone
– extra heavy duty version (XXL)*	200kg/31.5 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	100-240v (50-60Hz)
Maximum staircase angle	56° (52° with hinge)

† +45mm for extra heavy duty (XXL) 2 way swivel with extended arm spacers.

†† +20mm with arm spacers, +44mm with extended arm spacers.

^Not available for extra heavy duty (XXL). *Not available with hinge.



User support: 08700 118282

Handicare, 82 First Avenue, Pensnett Estate, Kingswinford, West Midlands DY6 7FJ, United Kingdom

Textphone: (180015) 020 7490 8913 ext 05 Voicephone: (18002) 020 7490 8913 ext 05

Global Dealer Support: +44 (0) 844 225 3121 enquiries@handicare.co.uk www.handicare.co.uk

Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



handicare



Handicare 1000

Straight stairlift



handicare

The benefits of a Handicare 1000 straight stairlift

If you are thinking about installing a stairlift the Handicare 1000 offers you a safe and reliable solution.

Handicare stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Handicare 1000 stairlift offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

With a choice of six upholstery colours the Handicare 1000 will compliment most colour schemes.

The slimline track is fitted to the stairs, not the wall, so it is unlikely you will need to redecorate following the installation of your Handicare 1000 stairlift.

Like all Handicare stairlifts the 1000 Series has a battery backup which means it will continue to work in the event of a power cut so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the slim track maximises clear space on the stairway and can be supplied with a manual or powered hinged track if required.

The Handicare 1000 also comes with two remote controls which means, once folded, you can easily send it to the upstairs charge point freeing up space in your hallway.



Safety as standard

Your Handicare 1000 stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Handicare 1000 stairlift has an easy release retractable seat belt for added peace of mind. If you prefer you can have an easy fasten Velcro seat belt.

Easy to use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Handicare 1000 will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.

