
**Handicare 2000 Simplicity**

<b>A</b> Standard swivel radius – with short arm option	720mm 685mm
<b>B</b> Minimum wall to seat back	96mm
<b>C</b> Seat depth	385mm
<b>D</b> Minimum folded width to footplate	400mm
<b>E</b> Minimum open width to edge of footplate	689mm
<b>F</b> Armrest width – external	635mm
<b>G</b> Armrest width – internal	480mm
<b>H</b> Seat back height	400mm
<b>I</b> Footplate to seat height	497-579mm*
<b>J</b> Minimum footplate height	70mm
Bottom parking space required at 45°	670mm
Minimum track intrusion into staircase	116mm
Standard track intrusion at bottom of stairs	500mm
Minimum track intrusion at bottom of stairs – reduced lower intrusion	420mm
Maximum weight limit – Heavy duty version	115kg/18 stone 137kg/21.5 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	100-240v (50-60Hz)
Maximum staircase angle	56°

\*Footplate to seat height reduced to 477mm if footplate is set to higher position (90mm).



BUTE

User support: 08700 118282

Handicare, 82 First Avenue, Pensnett Estate, Kingswinford, West Midlands DY6 7FJ, United Kingdom  
Textphone: (180015) 020 7490 8913 ext 05 Voicephone: (18002) 020 7490 8913 ext 05  
Global Dealer Support: +44 (0) 844 225 3121 enquiries@handicare.co.uk www.handicare.co.uk

Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



**Handicare 2000 Simplicity**  
Curved stairlift





# The benefits of the Simplicity seat and 2000 curved stairlift system from Handicare

If you are thinking about installing a stairlift, and your stairs turn corners or are spiral, Handicare's Simplicity seat and twin tube stairlift system offers you a safe and reliable solution.

Handicare's stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

When installed, the practical Simplicity seat it offers you all the standard functions that you would expect plus the choice to upgrade to a powered footplate so that folding the seat will automatically fold the footplate too. This means you do not need to bend to fold and unfold your footplate.

Supplied in a neutral shade (Bute) with wipe clean upholstery and a choice of three standard track colours (Light Grey, Gun Metal and Brown) the Simplicity seat and 2000 twin tube stairlift lift system will compliment most colour schemes.



All Handicare stairlifts have a battery backup which means that your stairlift will continue to work in the event of power cut so you do not need to worry about being stranded up or downstairs.

Handicare's twin tube stairlift system suits most staircases with bends or corners and offer's you the world's tightest internal bend. The compact unit folds up leaving the staircase free for other users and the bespoke track is designed to fit your individual staircase and maximise clear space on the stairway.

Fitted to the stairs, not the wall, it is unlikely you will need to redecorate following the installation of your Handicare stairlift. Your Handicare stairlift will also come with two remote controls which means that, once folded, you can easily send it to the upstairs charge point freeing up space in your hallway and on your stairs.



## Safety as standard

Your Handicare stairlift will come with safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs and will come with a retractable seatbelt as standard.

## Easy to use

The operating toggle is easy to use, simply hold it in the direction you wish to travel, up or down, and your Handicare stairlift will move smoothly and gently on its way.

