

Plinth 2000

Model 40 Neurology Couch

Overview:

Developed in conjunction with the Bobath Centre in London, our 'double' couch provides probably the most firm and stable treatment surface available. Being 102cm wide and lowering to just 42cm it is perfect for use in hospital gyms & physiotherapy departments. Now safely lifts 320kg from its lowest point.




Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Hydraulic or electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 320kg
- Retractable castors
- Hoist clearance 140mm
- Breathing hole

Overall Dimensions

Length:	189cm (74") approx
Width:	102cm (40") approx
Height:	42cm - 91cm approx (17" - 36") approx

Upholstery Dimensions



Recommended Options:



Foot Switch



Battery Back-up



Folding Cotside



Motor Back



Large Castor Base



Paper Roll Holder

NHS Approved Supplier



UK Made

t: +44 (0)1449 767887
e: sales@plinth2000.com
www.plinth2000.com

Operational Parameters



For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Model 40E - Electric version

Model 40H - Hydraulic version

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



White



Dove



Black



Chocolate



Damson



Wine



Red



Terracotta



Rose



Night



Ocean



Sky



Slate



Aqua



Green



Sage



Lime



Yellow



Beige

Plinth 2000 Ltd.
Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk, IP14 5PP

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com
www.plinth2000.com

Plinth | 2000