


Handicare: 130 years of lifting experience


With 700.000 satisfied customers worldwide and more than 130 years of experience, Handicare stands for high quality: quality in terms of safety, innovation, comfort and design. This is our guarantee that all of our products meet the highest possible standards in terms of safety and quality.

Handicare is proud to have the ISO 9001: 2015 quality certificate and in addition to this, all our products are tested and approved by worldwide independent inspection institutions such as TNO and Liftinstituut.

With offices around the globe, including the United Kingdom, Netherlands, Germany, France, Italy, Norway, China, United States and Canada, Handicare is a brand known for its quality, innovation and service.

Your local Handicare Business Partner



 Handicare Stairlifts B.V.
Newtonstraat 35
1704 SB Heerhugowaard
The Netherlands

 stairlifts@handicare.com

 www.handicare-stairlifts.com

Handicare 1100

Safe, Slim and Smart



NEW!



"Thanks to my stairlift,
I can live independently
in the home I love"



Safe: strong and reliable

The robust 1100 is powered by Handicare's patented four wheel drive technology with four quiet and reliable motors.

The lift is fully certified by the international stairlift safety standard EN 81-40 and designed for years and years of carefree daily use.

Thanks to the continuous charge functionality you can park the seat anywhere you like on the rail: the batteries are always charged and ready to go.



Slim: stylish compact design

The Handicare 1100 is ultra slim and compact and can be placed close to the wall, leaving a lot of free space on the staircase for other stair users to safely climb the stairs. It's the slimmest straight stairlift rail on the market!

To save even more space, the stairlift seat, footrest and armrests can be folded when not in use.

The seat will swivel at the top landing, making sure you can get on and off the stairlift easily and safely. Also, the 1100 comes as standard with safety sensors that will stop the lift in case of obstruction.



Handicare 1100: the slimmest straight stairlift

Developed with 130 years of experience, Handicare is proud to introduce the unique Handicare 1100 stairlift. The 1100's strong yet super slimline rail is the narrowest on the market so it obscures less of your stairs than any other straight stairlift system.

With its compact design and grease free rail, the 1100 blends in nicely with your interior and is a breeze to keep clean.

Smart: clean and tidy

Due to the patented technology, there is no tooth rack on the Handicare 1100. This means that unlike other stairlifts, there is no grease or oil on the rail collecting dust and dirt, making the rail easy to clean and safe to touch. Your staircase will stay spick and span!

